

Neobiotech EB System

Contributed by Neobiotech
Wednesday, 20 July 2011

Neobiotech EB System Â Neobiotech CMI implant EB external fixture

Â 1. bioseal headBone loss will be minimized through rounded switching form that makes enough tissue sealing in narrow space (0.5mm). 2. Taper-straight-taper bodyBy mixing straight and taper body`s merits, Sinus Quick has taper-straight-taper body type. Thus, drilling and implant ion are much easier. The body is designed to endure bite forces and tension. 3. Microgroove-MacrothreadThread has deep with shallowcurved shape to maximize the bone volume between threads and perform easy tapping. The body is designed to endure vertical bite force and lateral pressure effectively. 4. Magic ThreadThread has deep with shallow curved shape to maximize the bone volume between threads and perform easy tapping. The body is designed to endure vertical bite force and lateral pressure effectively. 5. SurfaceThe surface is treated by RBM(Resorbable Blast Media, Ra 1.2~1.5) with long stability. 6.Â ApexPowerful and well-defined threads exist at the peak of the apex. Drilling power is remarkable and initial fixation at the apex is excellent.